




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Water Bath with PC Tub, 11 l, 1200 W

€793.50
plus VAT &
Shipping

Product Images



Description

The inexpensive alternative!

Particularly inexpensive water baths with PC tank, heat resistant up to 120°C.

Temperature range:

5° above RT up to 100°C

Accuracy: approx. $\pm 0.3^\circ\text{C}$

Safety thermostat to switch off the unit in case of water shortage or overheating.

Connection: 220 Volt/50 HZ

CE-tested, easy handling, easy to clean.

standard scope of delivery:

- Suspended thermostat with connection cable and plug
- built-in safety thermostat
- seamless, transparent polycarbonate tray
- Bottom insert plate for direct setting of vessels or racks

Additional Information

No.	E-4090
Manufacturer (Brand)	neoLab
EAN	4058072005849
Transport temperature	Room temperature
Temp MAX	100 °C
Wide	270 mm
Height	160 mm
Length	180 mm
Depth	180 mm
Operation via	Knobs/Buttons
Voltage	230 V
Water baths and thermostats properties	with overheating protection cap/lid available as accessory
Application thermostat (temp. range)	Heat thermostat
Resolution °C	1 °C
Bath depth	160 mm
Bath opening	180 x 270 mm
Energy supply	Power grid
Heating capacity	1200 W
Material bath	Polycarbonate (PC)
Material inside	Polycarbonate (PC)
Temp MIN above ambient temperature	5 °C
Temp. accuracy	+/- 0.3 °C
Type thermostat (construction)	Immersion thermostat
Volume water bath	11 l

