



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Water Bath with PC Tub, 30.0 l, 2000 W

€1,062.00
plus VAT &
Shipping

Product Images



Description

The inexpensive alternative!

Particularly inexpensive water baths with PC tank, heat resistant up to 120°C.

Temperature range:

5° above RT up to 100°C

Accuracy: approx. $\pm 0.3^\circ\text{C}$

Safety thermostat to switch off the unit in case of water shortage or overheating.

Connection: 220 Volt/50 HZ

CE-tested, easy handling, easy to clean.

standard scope of delivery:

- Suspended thermostat with connection cable and plug
- built-in safety thermostat
- seamless, transparent polycarbonate tray
- Bottom insert plate for direct setting of vessels or racks

Additional Information

No.	E-4092
Manufacturer (Brand)	neoLab
EAN	4058072005863
Transport temperature	Room temperature
Temp MAX	100 °C
Wide	510 mm
Height	170 mm
Length	320 mm
Depth	320 mm
Operation via	Knobs/Buttons
Voltage	230 V
Water baths and thermostats properties	with overheating protection cap/lid available as accessory
Resolution °C	1 °C
Application thermostat (temp. range)	Heat thermostat
Energy supply	Power grid
Heating capacity	2000 W
Material bath	Polycarbonate (PC)
Bath depth	170 mm
Bath opening	320 x 510 mm
Temp. accuracy	+/- 0.3 °C
Volume water bath	30 l
Material inside	Polycarbonate (PC)
Temp MIN above ambient temperature	5 °C
Type thermostat (construction)	Immersion thermostat

