



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Corning water bath, digital; 6 liters, 230V, EU&UK power plug

€807.06
plus VAT &
Shipping

Product Images



Description

The personalized water bath made of corrosion-resistant stainless steel. The seamless tub holds up to a maximum of 6 liters, filling with dry bath balls is also possible. The temperature is indicated by the digital display, the temperature range is between RT+5 °C and 99°C. Heat loss is reduced by the included slanted lid. The water condensed on the lid does not drip down directly, but runs along the inside of the slanted lid.

Additional Information

No.	CQ-0228
Manufacturer (Brand)	Corning
old neoLab article no.	C-2050
EAN	4058072225773
Transport temperature	Room temperature
Temp MAX	99 °C
Wide	290 mm
Height	330 mm
Length	380 mm
Depth	380 mm
Operation via	Membrane keypad
Weight	8.6 kg
Voltage	230 V
Water baths and thermostats properties	with overheating protection with closure/lid
Resolution °C	0.1 °C
Application thermostat (temp. range)	Heat thermostat
Energy supply	Power grid
Heating capacity	1200 W
Material bath	Stainless steel
Bath depth	180 mm
Bath opening	170 x 320 mm
Temp. accuracy	+/- 0.1 °C
Volume water bath	6 l
Material inside	Stainless steel
Temp MIN above ambient temperature	5 °C
Type temperature control	PID
Type thermostat (construction)	Bath thermostat

