



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Julabo CORIO CD-200F Refrigerated Circulator

€3,730.00
plus VAT &
Shipping

Product Images



Description

Refrigerated circulators of the new CORIO series are characterized by a particularly good price/performance ratio. They are suitable for standard tasks and routine work in laboratories and industry. Bath vessel: stainless steel Bath cover: integrated Medium bath opening (cm/WxD/BD): 13 x 15 / 15 Pump function: pressure pump Pump type: submersible pump Product features: - Models for internal and external applications - Bright, white, highly visible display - Particularly quiet - Convenient pump switching between internal and external circulation - External pump connections - USB connection - Bath cover and drain cock are included - Removable ventilation grille - Chillers without side vents - Class III (FL) according to DIN 12876-1 - High-quality stainless steel bath vessels with bath cover and drain cock - More usable space in the bath due to space-saving evaporator coil Scope of delivery: 2 hose olives each suitable for hose 8 and 12 mm I.W. (pump connections with M16x1 external thread) Cooling capacity: Capacity data measured according to DIN 12876. Cooling capacities up to 20 °C measured with ethanol, above 20 °C measured with thermal oil, unless otherwise stated. Performance data apply at an ambient temperature of 20 °C. Performance values may deviate with other temperature control fluids.

Additional Information

No.	JN-0389
Manufacturer (Brand)	Julabo
VGKL number	103592701
EAN	4058072846893
Transport temperature	Room temperature
Net weight	26.30000 kg
Country of origin	Germany
Customs tariff number	84198998

