



Julabo DYNEO DD immersion thermostat

Product Images



Description

DYNEO DD immersion circulators stand for maximum flexibility. They are supplied with a stable and robust bath clamp and can be mounted on any bath vessels with up to 50 liters. The immersion circulators are used for internal and external temperature control and cover a working temperature range from +20 °C to +200 °C. The multilingual 3.5 inch color display and unique rotary knob provide easy and intuitive operation. Cooling coil: Accessories Pump function: Pressure pump Pump type: Submersible pump Product features: - For internal and external applications - Infinitely adjustable, powerful pressure pump - Flow rate 27 l/min, discharge pressure 0.7 bar - Easy switching between internal and external circulation - Large color TFT display, multilingual user interface - Easy operation via central rotary knob (controller) - Integrated programmer - External Pt100 connection integrated - USB connection - RS232 interface (optional) or analog interfaces (optional) - Flexible mounting - For internal and external applications Performance data: For temperature control applications near or below ambient temperature: use cooling coil or JULABO immersion cooler.

Additional Information

No.	JN-0658
Manufacturer (Brand)	Julabo
VGKL number	103590100
EAN	4058072847104
Net weight	2.50000 kg
Country of origin	Deutschland
Customs tariff number	84198998

