



Julabo DYNEO DD-BC26 circulation thermostat

Product Images



Description

DYNEO DD circulators for internal and external applications are equipped with closed bath vessels. The bath vessels are insulated to a high standard and have an integrated cooling coil for counter-cooling. The integrated drain cock allows safe and clean draining. The multilingual 3.5 inch color display and unique rotary knob provide easy and intuitive operation. Bath vessel: stainless steel Bath cover: integrated Medium bath opening (cm/WxD/BD): 26 x 35 / 20 Pump function: pressure pump Pump type: submersible pump Product features: - For internal and external applications - Infinitely adjustable, powerful pressure pump - Flow rate 27 l/min, delivery pressure 0.7 bar - Easy switching between internal and external circulation - Large color TFT display, multilingual user interface - Easy operation via central rotary knob (controller) - Integrated programmer - External Pt100 connection integrated - USB connection - RS232 interface (optional) or analog interfaces (optional) - Integrated drain cock for easy and safe draining - High-quality thermal insulation of the bath tank - For internal and external applications Performance data: For temperature control applications near or below ambient temperature: use cooling coil or JULABO immersion cooler. Scope of delivery: 2 hose olives each suitable for 8 and 12 mm I.w. tubing (pump connections with M16x1 male thread)

Additional Information

No.	JN-0670
Manufacturer (Brand)	Julabo
VGKL number	103591526
EAN	4058072847227
Net weight	18.70000 kg
Country of origin	Deutschland
Customs tariff number	84198998

