



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Immersion thermostat, stainless steel, 108 x 120 x 435 mm

€227.00
plus VAT &
Shipping

Product Images



Description

The shapely **immersion circulator** convinces with its compact shape and easy handling: It can be attached to containers with a maximum volume of 25 liters, can be programmed automatically, and is quickly cleaned after use and stored in a space-saving manner.

- User-friendly menu
- Heating time programming, automatic shut-off
- High accuracy, temperature stability: $\pm 0.1^{\circ}\text{C}$ (at 60°C)
- Alarm in case of overheating and low water level
- Protected design

Additional Information

No.	4-4423
Manufacturer (Brand)	neoLab
EAN	4058072157708
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	silver
Material	Stainless steel
Temp MAX	90°C
Temp MIN	40°C
Length	435 mm
Operation via	Knobs/Buttons
Weight	1.6 kg
Voltage	230 V
Application thermostat (temp. range)	Heat thermostat
Energy supply	Power grid
Heating capacity	1000 W
Temp. accuracy	$\pm 0.1^{\circ}\text{C}$
Type thermostat (construction)	Immersion thermostat

