



Julabo FK30-SL Calibration Thermostat

Product Images



Description

The calibration thermostats on this category are also suitable for calibration tasks down to -30 °C thanks to the additional integrated refrigeration machine. The bath cover with openings and Viton® feedthroughs is already included in the scope of delivery for all models. Cooling coil: not available Bath vessel: stainless steel Bath cover: integrated Average bath opening (cm/WxD/BT): 12 Pump type: submersible pump Product features: - LCD DIALOG DISPLAY for interactive operator guidance in plain text - ICC temperature control (Intelligent Cascade Control), self-optimizing - TCF Temperature Control features for optimization of the control behavior - Integrated programmer for 6 x 60 program steps - Highest temperature constancy up to ±0.005 °C (±0.02 °C at +200 °C) - Temperature homogeneity better than ±0.01 °C - Pt100 external sensor for reference temperature measurements optionally available - Temperature chamber with level-limited overflow Scope of delivery: Including bath cover with seven openings Cooling capacity: Performance data measured according to DIN 12876. Cooling capacities up to 20 °C measured with ethanol, above 20 °C measured with thermal oil, unless otherwise stated. Performance data apply at an ambient temperature of 20 °C. Performance values may deviate with other temperature control fluids.

Additional Information

No.	JN-0457
Manufacturer (Brand)	Julabo
EAN	4058072847937
Net weight	48.00000 kg
Temp MAX	200 °C
Temp MIN	-30 °C
Country of origin	Deutschland
Customs tariff number	84198998

