



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



COY O2 Control Glove Box, 1 person

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Our product specialist Juergen Traunecker is always available for further information.

Direct contact:

Tel.: 0174-3378337

Mail: j.traunecker@neolab.de

You want to work with your hypoxic cell cultures in a permanently controlled O₂ concentration?

The solution for this are the hypoxic glove boxes from Coylab.

Long-term incubation is performed in the chamber's humidity-regulated incubator.

O₂, CO₂ and temperature controllers allow direct monitoring of the chamber atmosphere.

Cell treatment in the largely sterile atmosphere is also possible without gloves (option).

Additional Information

No.	7-2653
Manufacturer (Brand)	COY
EAN	4058072194673
Transport temperature	Room temperature
Material	Polycarbonate (PC)
Temp MAX	45 °C
Temp MIN	3 °C
Wide	610 mm
Width/inside	590 mm
Height	720 mm
Height inside	610.0 mm
Length	1370 mm
Depth	1370 mm
Depth inside	1040 mm
Length inside	1040 mm
Voltage	230 V
Workbenches and work tents features	Underframe available as accessory
Worktop inside	0.61 m ²
Shelves inside	Shelf 4 tiers. variable
Number of workstations	1
Material chamber	Polycarbonate (PC)
Material inside	Polycarbonate (PC)
Port execution	without arm cuff
Type working chamber	Hypoxic working chamber
Type atmosphere control	N2 / CO2 / O2

