



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Centr.5430, Rotor A-2-MTP, knobs, 230V

**€4,776.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Centrifuge 5430, knobs, non-cooled, Low-Throughput Plate Solution, with rotor A-2-MTP, 230 V/50 - 60 Hz (EU) - Low-Throughput Plate Solution - knobs - with rotor A-2-MTP - non-cooled - 230 V, 50 - 60 Hz - 5 program keys for easy access to routine programs - Stores up to 50 user-defined settings - Automatic rotor and imbalance detection for maximum safety - Eppendorf QuickLock® system for fast opening and closing of rotor lid

## Additional Information

---

No.	VB-2808
Manufacturer (Brand)	Eppendorf
EAN	4058072672973
Gross weight	38.8
HSNumber	84211920
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	33.20000 kg
Speed MAX	17500 rpm
Volume MAX	300 mL
Country of origin	USA
Customs tariff number	84211920
Voltage	230 volt
Speed MIN	100 rpm
RZB	30130 g

