



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Centrifuge 5425, non-cooled, with rotor FA-24x2, membrane keyboard

Special Price

€2,288.80 was

€2,861.00

**plus VAT &
Shipping**

Product Images



Description

With its capacity for 24 tubes and a speed of up to $21,330 \times g$, the **Centrifuge 5425** is the new laboratory standard. As a device specifically for standard applications, it simply belongs in a modern laboratory, where it ensures optimal separation. The Centrifuge 5425 serves as a reliable partner, especially for high-throughput users who expect the highest ergonomic comfort, and offers particularly high flexibility with its six rotor options. It is perfectly suited for all modern molecular biology applications in Eppendorf and PCR tubes.

Additional Information

No.	VB-1828
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0811
EAN	4058072418809
Transport temperature	Room temperature
Speed MAX	15600 rpm
Color	White
sterile	No
Wide	390 mm
Height	240 mm
Length	240 mm
Depth	240 mm
Centrifuges properties	with timer with interface with continuous function including rotor(s) rotor is autoclavable rotor is exchangeable (Additional) rotors are available as accessories
Number of speed levels	> 10
Number of rotors supplied	1
Operation via	Membrane keypad
Speed MIN	2000 rpm
Fixed angle rotor	45°
Weight	15.6 kg
Power	280 W
Material rotor	Aluminum
Model centrifuge	Eppendorf. 5425
RZB	21130 g
Rotor(s) is/are suitable for	Reaction tubes. 24 x 1.5/ 2.0 ml
Voltage	230 V
Type rotor	Fixed angle rotor

Type centrifuge	Table centrifuge (small)
for ambient temperature	10 to 40 °C

