



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf (IVD) Centrifuge 5702 R, cooled, with rotor A-4-38, adapter for 15/50 mL conicals, rotary knobs

€6,147.00
plus VAT &
Shipping

Product Images



Description

The all-purpose **5702** low-speed **centrifuge family** is specifically designed for cell culture and clinical research laboratories with low to medium throughput. They fit easily under any bench thanks to their compact design and are quiet in operation, benefiting your working environment. They can be operated with six different rotor options to accommodate almost all vessel types.

' Cooled model **Centrifuge 5702 R** for heat-sensitive samples

Ergonomics

- ' Requires little floor space and fits on any lab bench
- ' Very quiet operation for a better working environment
- ' Low access height for easy loading and unloading of samples

Temperature management

- ' FastTemp function for fast pre-cooling (5702 R/RH)
- ' ECO shut-off switches on after 8 hours without operation to reduce energy consumption and extend compressor life (5702 R/RH)

Quality and safety

- ' Rotor chamber made of stainless steel, rust-proof and easy to clean
- ' Unbalance detection for maximum safety
- ' Key lock prevents unintentional adjustment of parameters
- ' Aerosol-tight caps for safe centrifugation of hazardous samples optionally available

Additional Information

No.	VB-2055
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-1819
EAN	4058072420925
Transport temperature	Room temperature
Speed MAX	4400 rpm
Color	White
Temp MAX	42 °C
Temp MIN	-9 °C
sterile	No
Wide	581 mm
Height	270 mm
Length	381 mm
Depth	381 mm
Operation via	Knobs/Buttons
Weight	35.1 kg
Voltage	230 V
for ambient temperature	10 to 40 °C
Centrifuges properties	programmable with timer with cooling with continuous function including rotor(s) rotor is autoclavable rotor is exchangeable (Additional) rotors are available as accessories
Speed MIN	100 rpm
Number of speed levels	> 10
Power	380 W
Number of rotors supplied	1
Material rotor	Aluminum
Maximum relative humidity at 20 °C	75.0 %
Model centrifuge	Eppendorf. 5702 R

RZB	3000 g
Type rotor	Swing-out rotor
Type centrifuge	Table centrifuge (large)

