



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf (IVD) benchtop centrifuge 5804, without rotor

€4,500.00
plus VAT &
Shipping

Product Images



Description

The **5804 centrifuge** is a high-speed centrifuge for medium throughput. It is suitable for molecular biology applications with tubes up to 250 mL and offers more flexibility with additional swing-out and fixed-angle rotors as well as the inclusion of deepwell plates.

Additional Information

No.	VB-2061
Manufacturer (Brand)	Eppendorf
old neolab article no.	E-1892
EAN	4058072420987
Transport temperature	Room temperature
Speed MAX	14000 rpm
Color	White
Wide	550 mm
Height	340 mm
Length	470 mm
Depth	470 mm
Operation via	Membrane keypad
Weight	55 kg
Voltage	230 V
for ambient temperature	2 to 40 °C
Centrifuges properties	programmable with timer with continuous function rotor is replaceable (additional) rotors are available as accessories
Speed MIN	200 rpm
Number of speed levels	> 10
Power	900 W
Maximum relative humidity at 20 °C	75.0 %
Model centrifuge	Eppendorf. 5804
RZB	20193 g
Type centrifuge	Table centrifuge (large)

