



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Centrifuge 5810 R, cooled, without rotor, membrane keyboard

€9,268.97
plus VAT &
Shipping

Product Images



Description

With its renowned quality and reliability, the Eppendorf **Centrifuge 5810/5810 R** is the most cost-effective solution for medium to high throughput applications - now and in the future. It combines exceptional flexibility and capacity for both tubes and plates with a particularly small footprint. With the new generation of rotors, the capacity of the 5810/5810 R centrifuge increases to max. 4 × 750 mL (or max. 28 × 50 mL/56 × 15 mL). It is currently the most compact 3-liter centrifuge on the market.

Fits in any laboratory

- ' The most compact 3 liter centrifuge on the market
- ' Can be used under shelves thanks to low lid opening height
- ' Quiet operation for a better working environment
- ' Low access height for easy loading and unloading
- ' Centrifuge lid with soft-touch lid lock

Eppendorf ^{QuickLock®} caps for one-hand operation

- ' Significantly improved handling, easy one-handed operation
- ' Thanks to the handle, the closed carrier can be carried to the microbiological safety cabinet and opened there in a safe environment.
- ' Aerosol tightness has been verified by Public Health England, Porton Down, UK

Additional Information

No.	VB-2186
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0842
EAN	4058072422233
Transport temperature	Room temperature
autoclavable at 121°C	No
Speed MAX	14000 rpm
Color	White
Temp MAX	40 °C
Temp MIN	-9 °C
sterile	No
Wide	610 mm
Height	350 mm
Length	700 mm
Depth	700 mm
Operation via	Membrane keypad
Weight	99 kg
Voltage	230 V
Centrifuges properties	programmable with timer with cooling with continuous function rotor is replaceable (additional) rotors are available as accessories
Number of speed levels	> 10
Speed MIN	200 rpm
Power	1650 W
Model centrifuge	Eppendorf. 5810 R
RZB	20913 g
Type centrifuge	Table centrifuge (large)

