



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Centrifuge 5810 R, cooled, with rotor S-4-104, adapter for 13/16 mm round-bottom tubes, membrane keypad

**€12,246.70**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

With its renowned quality and reliability, the Eppendorf **Centrifuge 5810/5810 R** is the most cost-effective solution for medium to high throughput applications - now and in the future. It combines exceptional flexibility and capacity for both tubes and plates with a particularly small footprint. With the new generation of rotors, the capacity of the 5810/5810 R centrifuge increases to max. 4 × 750 mL (or max. 28 × 50 mL/56 × 15 mL). It is currently the most compact 3-liter centrifuge on the market.

### Fits in any laboratory

- ' The most compact 3 liter centrifuge on the market
- ' Can be used under shelves thanks to low lid opening height
- ' Quiet operation for a better working environment
- ' Low access height for easy loading and unloading
- ' Centrifuge lid with soft-touch lid lock

### Eppendorf <sup>QuickLock®</sup> caps for one-hand operation

- ' Significantly improved handling, easy one-handed operation
- ' Thanks to the handle, the closed carrier can be carried to the microbiological safety cabinet and opened there in a safe environment.
- ' Aerosol tightness has been verified by Public Health England, Porton Down, UK

## Additional Information

No.	VB-2193
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0846
EAN	4058072422301
Transport temperature	Room temperature
autoclavable at 121°C	No
Speed MAX	14000 rpm
Color	White
Temp MAX	40 °C
Temp MIN	-9 °C
sterile	No
Wide	610 mm
Height	350 mm
Length	700 mm
Depth	700 mm
Operation via	Membrane keypad
Weight	99 kg
Voltage	230 V
Centrifuges properties	programmable with timer with cooling with continuous function including rotor(s) rotor is autoclavable rotor is exchangeable (Additional) rotors are available as accessories
Speed MIN	200 rpm
Number of speed levels	> 10
Power	1650 W
Number of rotors supplied	1
Material rotor	Aluminum
Model centrifuge	Eppendorf. 5810 R
RZB	20913 g

Type rotor

Swing-out rotor

Type centrifuge

Table centrifuge (large)

