



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Eppendorf (IVD) Centrifuge 5810, non-cooled,
with rotor A-4-81, adapter for 15/50 mL conicals,
membrane keypad**

**€8,860.00
plus VAT &
Shipping**

Product Images



Description

With its renowned quality and reliability, the Eppendorf **Centrifuge 5810/5810 R** is the most cost-effective solution for medium to high throughput applications - now and in the future. It combines exceptional flexibility and capacity for both tubes and plates with a particularly small footprint. With the new generation of rotors, the capacity of the 5810/5810 R centrifuge increases to max. 4 × 750 mL (or max. 28 × 50 mL/56 × 15 mL). It is currently the most compact 3-liter centrifuge on the market.

Fits in any laboratory

- ' The most compact 3 liter centrifuge on the market
- ' Can be used under shelves thanks to low lid opening height
- ' Quiet operation for a better working environment
- ' Low access height for easy loading and unloading
- ' Centrifuge lid with soft-touch lid lock

Eppendorf ^{QuickLock®} caps for one-hand operation

- ' Significantly improved handling, easy one-handed operation
- ' Thanks to the handle, the closed carrier can be carried to the microbiological safety cabinet and opened there in a safe environment.
- ' Aerosol tightness has been verified by Public Health England, Porton Down, UK

Additional Information

No.	VB-2135
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-1994
EAN	4058072421724
Transport temperature	Room temperature
autoclavable at 121°C	No
Speed MAX	14000 rpm
Color	White
sterile	No
Wide	610 mm
Height	350 mm
Length	540 mm
Depth	540 mm
Operation via	Membrane keypad
Weight	68 kg
Voltage	230 V
Centrifuges properties	programmable with timer with continuous function including rotor(s) rotor is autoclavable rotor is replaceable (Additional) rotors are available as accessories
Speed MIN	200 rpm
Number of speed levels	> 10
Power	900 W
Number of rotors supplied	1
Material rotor	Aluminum
Model centrifuge	Eppendorf. 5810
RZB	20913 g
Type rotor	Swing-out rotor
Type centrifuge	Table centrifuge (large)

