



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf (IVD) Centrifuge 5810 R, cooled, with rotor S-4-104, adapter for 15/50 mL conicals, membrane keypad

**€12,854.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

With its renowned quality and reliability, the Eppendorf **Centrifuge 5810/5810 R** is the most cost-effective solution for medium to high throughput applications - now and in the future. It combines exceptional flexibility and capacity for both tubes and plates with a particularly small footprint. With the new generation of rotors, the capacity of the 5810/5810 R centrifuge increases to max. 4 × 750 mL (or max. 28 × 50 mL/56 × 15 mL). It is currently the most compact 3-liter centrifuge on the market.

### Fits in any laboratory

- ' The most compact 3 liter centrifuge on the market
- ' Can be used under shelves thanks to low lid opening height
- ' Quiet operation for a better working environment
- ' Low access height for easy loading and unloading
- ' Centrifuge lid with soft-touch lid lock

### Eppendorf <sup>QuickLock®</sup> caps for one-hand operation

- ' Significantly improved handling, easy one-handed operation
- ' Thanks to the handle, the closed carrier can be carried to the microbiological safety cabinet and opened there in a safe environment.
- ' Aerosol tightness has been verified by Public Health England, Porton Down, UK

## Additional Information

No.	VB-2192
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-2003
EAN	4058072422295
Transport temperature	Room temperature
autoclavable at 121°C	No
Speed MAX	14000 rpm
Color	White
Temp MAX	40 °C
Temp MIN	-9 °C
sterile	No
Wide	610 mm
Height	350 mm
Length	700 mm
Depth	700 mm
Operation via	Membrane keypad
Weight	99 kg
Voltage	230 V
Centrifuges properties	programmable with timer with cooling with continuous function including rotor(s) rotor is autoclavable rotor is exchangeable (Additional) rotors are available as accessories
Number of speed levels	> 10
Number of rotors supplied	1
Speed MIN	200 rpm
Power	1650 W
Material rotor	Aluminum
Model centrifuge	Eppendorf. 5810 R
RZB	20913 g

Type rotor

Swing-out rotor

Type centrifuge

Table centrifuge (large)

