



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Centrifuge 5910 R, cooled, without rotor, membrane keyboard

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

The **5910 R centrifuge** combines a high capacity of up to 4 × 750 mL with exceptional flexibility. Its main swing-out rotor accepts both conical tubes and plates – no need to exchange rotor buckets or adapters. This improves handling and saves time. It is equipped with a powerful cooling system with state-of-the-art temperature management so you can safely centrifuge your samples.

### **Low access height for easy operation**

The 5910 R centrifuge has been designed with a particularly low access height. This not only allows easy loading and unloading, but everyone in the lab can also easily reach the open centrifuge lid.

### **Softtouch lid closure**

No unnecessary force is required when closing the centrifuge lid. The soft-touch lid closure makes it easy to close the centrifuge.

### **State-of-the-art operating system and temperature management**

' Menu-guided multilingual operation with backlight

' FastTemp<sup>pro</sup>® function for pre-cooling programmable by day and time; in combination with the ECO shut-off function, this significantly reduces power consumption

' Dynamic compressor control for precise cooling performance

## Additional Information

No.	VB-2366
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0847
EAN	4058072423995
Transport temperature	Room temperature
Speed MAX	14000 rpm
Color	White
Temp MAX	40 °C
Temp MIN	-11 °C
sterile	No
Wide	660 mm
Height	370 mm
Length	720 mm
Depth	720 mm
Operation via	Membrane keypad
Weight	109 kg
Voltage	230 V
Centrifuges properties	programmable with timer with cooling with continuous function rotor is replaceable (additional) rotors are available as accessories
Speed MIN	100 rpm
Number of speed levels	> 10
Power	1650 W
Maximum relative humidity at 20 °C	75.0 %
Model centrifuge	Eppendorf. 5910 R
RZB	22132 g
Type centrifuge	Table centrifuge (large)

