



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf (IVD) Centrifuge 5910 R, cooled, with rotor S-4x500, incl. rectangular beaker and adapter for 15/50 mL conicals, membrane keypad

**Price on
request**

**plus VAT &
Shipping**

Product Images



Description

The **5910 R centrifuge** combines a high capacity of up to 4 × 750 mL with exceptional flexibility. Its main swing-out rotor accepts both conical tubes and plates – no need to exchange rotor buckets or adapters. This improves handling and saves time. It is equipped with a powerful cooling system with state-of-the-art temperature management so you can safely centrifuge your samples.

Low access height for easy operation

The 5910 R centrifuge has been designed with a particularly low access height. This not only allows easy loading and unloading, everyone in the lab can also easily reach the open centrifuge lid.

Softtouch lid closure

No unnecessary force is required when closing the centrifuge lid. The soft-touch lid closure makes it easy to close the centrifuge.

State-of-the-art operating system and temperature management

' Menu-guided multilingual operation with backlighting

' FastTemp^{pro}® function for pre-cooling programmable by day and time; in combination with the ECO shut-off function, this significantly reduces power consumption

' Dynamic compressor control for precise cooling performance

Swing-out rotor with exceptional versatility

With the large S-4xUniversal swing-bucket rotor, tubes and plates can be used simultaneously by replacing the adapter instead of the rotor bucket.

Your advantages:

' Time-saving: no need to exchange buckets

' Cost-saving: no need for separate plate buckets

' Space-saving: no need to store different buckets

In addition, selected conical tubes and plate formats can be used with the three new innovative 'universal adapters':

' Eppendorf^{Tubes}® 5.0 mL/15 mL conical/DWP

' 50 mL conical/MTP

' 50 mL conical/250 mL/MTP

Your convenience:

Even more convenience when using conical tubes and plates.

Additional Information

| | |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| No. | VB-2376 |
| Manufacturer (Brand) | Eppendorf |
| old neoLab article no. | E-2011 |
| EAN | 4058072424091 |
| Transport temperature | Room temperature |
| Speed MAX | 14000 rpm |
| Color | White |
| Temp MAX | 40 °C |
| Temp MIN | -11 °C |
| sterile | No |
| Wide | 660 mm |
| Height | 370 mm |
| Length | 720 mm |
| Depth | 720 mm |
| Operation via | Membrane keypad |
| Weight | 109 kg |
| Voltage | 230 V |
| Centrifuges properties | programmable with timer with cooling with continuous function including rotor(s) rotor is autoclavable rotor is exchangeable (Additional) rotors are available as accessories |
| Speed MIN | 100 rpm |
| Number of speed levels | > 10 |
| Power | 1650 W |
| Number of rotors supplied | 1 |
| Material rotor | Aluminum |
| Maximum relative humidity at 20 °C | 75.0 % |
| Model centrifuge | Eppendorf. 5910 R |
| RZB | 22132 g |

Type rotor

Swing-out rotor

Type centrifuge

Table centrifuge (large)

