



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Centrifuge 5920 R, cooled, without rotor, membrane keyboard

**€12,420.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The **5920 R centrifuge** offers an exceptionally large capacity in a compact and ergonomic product design. It has similar dimensions to refrigerated competitor models in the 3-liter segment, but has an outstanding capacity of up to 4 × 1,000 mL. This makes it the ideal unit for high-throughput applications. It is equipped with a powerful cooling system with state-of-the-art temperature management, so you can centrifuge your samples safely.

### Modern operating system

The Centrifuge 5920 R features a modern operating system with excellent functionality

- ' Menu-driven, multilingual operating menu and large, backlit display
- ' Stores up to 99 user-defined programs
- ' 5 program keys for quick access to routine programs
- ' FastTemp pro<sup>®</sup> function for automatic pre-cooling at preprogrammable time

## Additional Information

No.	VB-2560
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0854
EAN	4058072425692
Transport temperature	Room temperature
Speed MAX	13700 rpm
Color	White
Temp MAX	40 °C
Temp MIN	-11 °C
sterile	No
Wide	710 mm
Height	400 mm
Length	740 mm
Depth	740 mm
Operation via	Membrane keypad
Weight	139 kg
Voltage	230 V
Centrifuges properties	programmable with timer with cooling with continuous function rotor is replaceable (additional) rotors are available as accessories
Speed MIN	100 rpm
Number of speed levels	> 10
Power	1650 W
Maximum relative humidity at 20 °C	75.0 %
Model centrifuge	Eppendorf. 5920 R
RZB	21194 g
Type centrifuge	Table centrifuge (large)

