



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf (IVD) Centrifuge 5920 R, cooled, without rotor, membrane keyboard

€13,035.60
plus VAT &
Shipping

Product Images



Description

The **5920 R centrifuge** offers an exceptionally large capacity in a compact and ergonomic product design. It has similar dimensions to refrigerated competitor models in the 3-liter segment, but has an outstanding capacity of up to 4 × 1,000 mL. This makes it the ideal unit for high-throughput applications. It is equipped with a powerful cooling system with state-of-the-art temperature management, so you can centrifuge your samples safely.

Modern operating system

The 5920 R centrifuge features a modern operating system with excellent functionality

- ' Menu-driven, multilingual operating menu and large, backlit display
- ' Stores up to 99 user-defined programs
- ' 5 program keys for quick access to routine programs
- ' FastTemp pro[®] function for automatic pre-cooling at preprogrammable time

Additional Information

No.	VB-2561
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-2013
EAN	4058072425708
Transport temperature	Room temperature
Speed MAX	13700 rpm
Color	White
Temp MAX	40 °C
Temp MIN	-11 °C
sterile	No
Wide	710 mm
Height	400 mm
Length	740 mm
Depth	740 mm
Operation via	Membrane keypad
Weight	139 kg
Voltage	230 V
Centrifuges properties	programmable with timer with cooling with continuous function rotor is replaceable (additional) rotors are available as accessories
Speed MIN	100 rpm
Number of speed levels	> 10
Power	1650 W
Maximum relative humidity at 20 °C	75.0 %
Model centrifuge	Eppendorf. 5920 R
RZB	21194 g
Type centrifuge	Table centrifuge (large)

