



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Eppendorf Centr.5920R, rotor S-4x750, adapt.
15/50mL, with cooling, membrane keyboard,
230V**

**€15,390.00
plus VAT &
Shipping**

Product Images



Description

Centrifuge 5920 R G, 230 V, 50-60 Hz, incl. rotor S-4x750, round beaker and adapter 15 mL/50 mL conical tubes

Additional Information

No.	VB-2568
Manufacturer (Brand)	Eppendorf
EAN	4058072425777
Gross weight	169.2
HSNumber	84211920
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	152.00000 kg
Speed MAX	13700 rpm
Country of origin	Germany
Customs tariff number	84211920
Speed MIN	100 rpm
RZB	21194 g

