



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Thermo Scientific Centrifuge Pico 17 + Rotor 24 x 1.5 ml, w. ClickSeal Biocontainment Lid

**€2,410.00**  
**plus VAT &**  
**Shipping**



## Description

---

Microcentrifuge Heraeus Pico 17, ventilated  
incl. rotor 3-1885

- generous capacity allows for a larger number of samples per run
- aerosol-tight ClickSeal™ rotor provides first-class safety and comfort
- fast acceleration and deceleration
- user-friendly controls
- high-quality materials allow vigorous cleaning and autoclaving

### Technical data

max. RZB:

17,000 x g

**rpm:**

13,300

noise level:

< 55 dBA

timer : 1- 99 min:1 min

dimensions H x W x D (mm) :

Pico17: 225 x 243 x 352

### **Order No. 3-1885 Rotor**

Rotor for 24 x 1.5/2.0 ml, with aerosol-tight ClickSeal biocontainment lid

## Additional Information

No.	LD-1919
Manufacturer (Brand)	Thermo Scientific
old neolab article no.	3-1860
EAN	4058072287665
Transport temperature	Room temperature
Speed MAX	13300 rpm
Color	White
Wide	352 mm
Height	225 mm
Length	243 mm
Depth	243 mm
Operation via	Membrane keypad
Weight	10.5 kg
Voltage	230 V
for ambient temperature	5 to 40 °C
Centrifuges properties	programmable with timer with continuous function including rotor(s) rotor is autoclavable rotor is replaceable (Additional) rotors are available as accessories
Speed MIN	800 rpm
Number of speed levels	> 10
Power	180 W
Number of rotors supplied	1
Rotor(s) is/are suitable for	Reaction tubes. 24 x 1.5/ 2.0 ml
Spare part	LD-1998
Fixed angle rotor	45°
Material rotor	Aluminum
Maximum relative humidity at 20 °C	80.0 %
Model centrifuge	Thermo Scientific Heraeus Pico 17
RZB	17000 g

Type rotor

Fixed angle rotor

Type centrifuge

Table centrifuge (small)

