



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Thermo Scientific Centrifuge Fresco 17 + Rotor 24 x 1.5 ml, w. ClickSeal Biocontainment Lid

**€5,800.00**  
**plus VAT &**  
**Shipping**

### Product Images

## Description

---

Microcentrifuge Heraeus Fresco 17, cooled  
incl. rotor 3-1885

- generous capacity allows for a larger number of samples per run
- aerosol-tight ClickSeal™ rotor provides first-class safety and comfort
- fast acceleration and deceleration
- user-friendly controls
- high-quality materials allow vigorous cleaning and autoclaving

### Technical data

max. RZB:

17,000 x g

**rpm:**

13.300

Noise level:

< 55 dBA

Timer: 1- 99 min:1 min

Temperature: -9 °C to +40 °C:1°C

Dimensions H x W x D (mm) :

Fresco 17: 330 x 295 x 445

### **Order No. 3-1885 Rotor**

Rotor for 24 x 1.5/2.0 ml, with aerosol-tight ClickSeal biocontainment lid

## Additional Information

No.	LD-1921
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1861
EAN	4058072287689
Transport temperature	Room temperature
Speed MAX	13300 rpm
Color	White
Temp MAX	40 °C
Temp MIN	-9 °C
Wide	445 mm
Height	330 mm
Length	295 mm
Depth	295 mm
Operation via	Membrane keypad
Weight	28 kg
Voltage	230 V
for ambient temperature	5 to 40 °C
Centrifuges properties	programmable with timer with cooling with continuous function including rotor(s) rotor is autoclavable rotor is exchangeable (Additional) rotors are available as accessories
Speed MIN	2000 rpm
Number of speed levels	> 10
Power	330 W
Number of rotors supplied	1
Rotor(s) is/are suitable for	Reaction tubes. 24 x 1.5/ 2.0 ml
Spare part	LD-1998
Fixed angle rotor	45°
Material rotor	Aluminum

Maximum relative humidity at 20 °C	80.0 %
Model centrifuge	Thermo Scientific Heraeus Fresco 17
RZB	17000 g
Type rotor	Fixed angle rotor
Type centrifuge	Table centrifuge (small)

