



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Centrifuge Fresco 21 + Rotor 24 x 1.5 ml, w. ClickSeal Biocontainment Lid

€6,490.00
plus VAT &
Shipping

Product Images



Description

Microcentrifuges Heraeus Fresco 21, cooled
incl. rotor

- generous capacity allows for a larger number of samples per run
- aerosol-tight ClickSeal™ rotor offers first-class safety and convenience
- double-row rotor eliminates the need for adapters
- fast acceleration and deceleration
- user-friendly controls
- high-quality materials allow vigorous cleaning and autoclaving

Technical data

max. RZB:

21,100 x g

rpm:

14.800

Noise level:

< 56 dBA

Timer: 1- 99 min:1 min

Temperature: -9 °C to +40 °C:1°C

Dimensions H x W x D (mm) :

Fresco 21: 330 x 295 x 445

Order No. 3-1863

Fresco 21 incl. rotor for 24 x 1.5/2.0 ml, with aerosol-tight ClickSeal biocontainment lid

Additional Information

| | |
|------------------------------------|---|
| No. | LD-1922 |
| Manufacturer (Brand) | Thermo Scientific |
| old neoLab article no. | 3-1863 |
| EAN | 4058072287696 |
| Transport temperature | Room temperature |
| Speed MAX | 14800 rpm |
| Color | White |
| Temp MAX | 40 °C |
| Temp MIN | -9 °C |
| Wide | 445 mm |
| Height | 330 mm |
| Length | 295 mm |
| Depth | 295 mm |
| Centrifuges properties | programmable with timer with cooling with continuous function including rotor(s) rotor is autoclavable rotor is exchangeable (Additional) rotors are available as accessories |
| Number of speed levels | > 10 |
| Number of rotors supplied | 1 |
| Operation via | Membrane keypad |
| Speed MIN | 2000 rpm |
| Spare part | LD-1998 |
| Fixed angle rotor | 45° |
| Weight | 28 kg |
| Power | 370 W |
| Material rotor | Aluminum |
| Maximum relative humidity at 20 °C | 80.0 % |
| Model centrifuge | Thermo Scientific Heraeus Fresco 21 |
| RZB | 21100 g |

| | |
|------------------------------|----------------------------------|
| Rotor(s) is/are suitable for | Reaction tubes. 24 x 1.5/ 2.0 ml |
| Voltage | 230 V |
| Type rotor | Fixed angle rotor |
| Type centrifuge | Table centrifuge (small) |
| for ambient temperature | 5 to 40 °C |

