



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand Bottle Attachment Burette Titrette 0-25 ml:0.001 ml

€924.50
plus VAT &
Shipping

Product Images



Description

With the bottle-top burette

Titrette® you titrate quickly and reliably with maximum precision.

High accuracy: Precise titration within the error limits of Class A for glass burettes.

Smooth-running gear: For drop-by-drop titration.

Compact design: Lightweight and stable.

Easy handling: No need to switch between "filling" and "titrating"

Easy maintenance: Cleaning and parts replacement possible in the laboratory.

Temperature range: +15°C to +40°C

Vapor pressure: up to 500 mbar

Viscosity: up to 500 mm²/s

Rel. humidity: 20 to 90 %

Accuracy:

Order No. 3-0816 : $\leq \pm 0.07$ %

Coefficient of variation:

Order No. 3-0816: ≤ 0.025 %

(Adjusted to "Ex". Based on the nominal volume of the instrument)

The SafetyPrime™ back-dosing valve enables fast, bubble-free venting with virtually no loss of media.

Scope of delivery

Adapter for GL 32, S 40 and NS 29/32

Titrette®, conformity certified, quality certificate, telescopic aspirating tube (170 to 330 mm), back dosing tube, 2 micro batteries, 2 brown colored light protection windows.

Optional: Mir interface RS 232

Accessories on request.

Additional Information

No.	BN-1441
Manufacturer (Brand)	Brand
old neolab article no.	3-0816
EAN	4058072208110
Transport temperature	Room temperature
Volume MAX	25 ml
Volume MIN	0 ml
Display (reading, setting)	digital
Resolution Bottle Top Dispenser	0 to 20ml: 0.001 ml 20 to 25ml: 0.01 ml
Functionality	electronic
Model Bottle Top Dispenser	Brand Titrette
Correctness at Vol. MAX	0.0007
Type Bottle Top Dispenser	Burette (digital)
Type certificate	Declaration of conformity (standardization) calibration certificate from manufacturer
Coefficient of variation at Vol. MAX	0.0003
for thread/DM	GL 32 GL 45 S 40
adjusted to	"ex"

