



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand Bottle Attachment Burette Titrette 0-50 ml:0.002 ml

€924.50
plus VAT &
Shipping

Product Images



Description

With the bottle-top burette

Titrette® you titrate quickly and reliably with maximum precision.

High accuracy: Precise titration within the error limits of Class A for glass burettes.

Smooth-running gear: For drop-by-drop titration.

Compact design: Lightweight and stable.

Easy handling: No need to switch between "filling" and "titrating"

Easy maintenance: Cleaning and parts replacement possible in the laboratory.

Temperature range: +15°C to +40°C

Vapor pressure: up to 500 mbar

Viscosity: up to 500 mm²/s

Rel. humidity: 20 to 90 %

Accuracy:

OrderNo. 3-0817: $\leq \pm 0.06$ %

Coefficient of variation:

Order No. 3-0817: ≤ 0.02 %

(Adjusted to "Ex". Based on the nominal volume of the device)

The SafetyPrime™ back-dosing valve enables fast, bubble-free venting with virtually no loss of media.

Scope of delivery

Adapter for GL 32, S 40 and NS 29/32

Titrette®, conformity certified, quality certificate, telescopic aspirating tube (170 to 330 mm), back dosing tube, 2 micro batteries, 2 brown colored light protection windows.

Optional: Mir interface RS 232

Accessories on request.

Additional Information

No.	BN-1442
Manufacturer (Brand)	Brand
old neoLab article no.	3-0817
EAN	4058072208127
Transport temperature	Room temperature
Volume MAX	50 ml
Volume MIN	0 ml
Display (reading, setting)	digital
Functionality	electronic
Model Bottle Top Dispenser	Brand Titrette
Resolution Bottle Top Dispenser	0 to 20 ml: 0.002 ml 20 to 50 ml: 0.01 ml
Correctness at Vol. MAX	0.0006
Type Bottle Top Dispenser	Burette (digital)
Type certificate	Declaration of conformity (standardization) calibration certificate from manufacturer
Coefficient of variation at Vol. MAX	0.0002
for thread/DM	GL 32 GL 45 S 40
adjusted to	"ex"

