



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® neoBuret e Digital burette with magnetic stirrer

€1,286.00
plus VAT &
Shipping

Product Images



Description

The **neoBuret e** is the first **digital burette** in the neoLabLine family. Easily and conveniently titrate the desired solution drop by drop to the desired result.

A high degree of accuracy and reproducibility is especially important in quantitative analysis. The **neoBuret e** is intuitive to operate and also delivers precise values for a satisfactory result. Acid-base reactions, precipitation titrations, redox reactions or even pharmacological dose adjustments have never been so easy.

- High-quality motor with excellent precision and accuracy
- Electronic control reduces the possibility of errors
- Remote control eliminates manual interference during operation
- With built-in magnetic stirrer

Additional Information

No.	7-1770
Manufacturer (Brand)	neoLab
EAN	4058072158491
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Acrylonitrile butadiene styrene (ABS)
Volume MAX	99.99 ml
Volume MIN	0.01 ml
Display (reading, setting)	graphic display
Functionality	electronic
Number of parts in set	3
Correctness at Vol. MAX	0.002
Type Bottle Top Dispenser	Burette (digital)
Type certificate	Declaration of conformity (standardization)
Coefficient of variation at Vol. MAX	0.0007
for thread/DM	GL 45

