



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hirschmann® opus® titration, 50 ml, electronic motorized bottle-top titrator with touchscreen module, 100-240 volts

€1,924.00
plus VAT &
Shipping

Product Images



Description

opus® titration, 50 ml, 100-240 Volt, (EU, US, UK plug), electronic motorized bottle top titrator with touch screen module "titration", with reagent feedback system, with declaration of conformity according to DIN 12600, serial or USB interface, complete with screw thread A 45 mm, 3 thread adapters (A 32, A 38 and S 40), 1 suction hose with union nut, 1 titration ejector unit and 1 power supply unit, Data Power cable 0.75 m, operating instructions, 1 valve mounting cap, with individual serial number and calibration certificate.

Additional Information

No.	NH-1628
Reg. trademark	opus ®
Manufacturer (Brand)	Hirschmann
VGKL number	103782050
EAN	4058072812058
Manufacturer EAN	4049338020485
Gross weight	1.845
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	1.16400 kg
Country of origin	Germany
Packing width	0.24500 m
Packing height	0.10800 m
Packing depth	0.36800 m
Packaging volume	0.00973728 cc
Customs tariff number	84798997

