



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



tec-lab temperature logger for CXR 500

€1,960.00
plus VAT &
Shipping

Product Images



Description

Temperature logger for cryo sand container LN2 Order No. C-7557 Suitable for CXR 500

Additional Information

No.	TL-0034
Manufacturer (Brand)	tec-lab
old neoLab article no.	C-7557
EAN	4058072373528
Transport temperature	Room temperature
Temp MAX	40 °C
Temp MIN	-199 °C
suitable for	Cryo sand container (LN2)
Type interface	USB
Figure	2
Measuring accuracy °C	3 °C

