



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## tec-lab temperature logger for CXR 100

**€1,358.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Temperature logger for cryo sand container LN2 Order No. C-7556 Suitable for CXR 100

## Additional Information

---

|                        |                           |
|------------------------|---------------------------|
| No.                    | TL-0035                   |
| Manufacturer (Brand)   | tec-lab                   |
| old neoLab article no. | C-7556                    |
| EAN                    | 4058072373535             |
| Transport temperature  | Room temperature          |
| Temp MAX               | 40 °C                     |
| Temp MIN               | -199 °C                   |
| suitable for           | Cryo sand container (LN2) |
| Type interface         | USB                       |
| Figure                 | 6                         |
| Measuring accuracy °C  | 3 °C                      |

