



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Amarell USB data logger, type 98581, for temperature with LCD display, measuring range

**Price on
request
plus VAT &
Shipping**



Description

USB data logger, type 98581, for temperature with LCD display, measuring range -40...+85 in 0,1°C, accuracy +/-0,6°C from -20 to +50°C, otherwise +/-1,2°C, 32000 memory values, measuring cycle 1sec. - 6h, complete with Windows software and battery, dimensions 136x30mm

Additional Information

No.	LL-0687
Manufacturer (Brand)	Amarell
VGKL number	101018581
EAN	4058072462925
Temp MAX	85 °C
Temp MIN	-40 °C
Country of origin	Taiwan, Province of China
Customs tariff number	90258040
measured size(s)	Temperature

