



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Coating Thickness Gauge TB 1000-0.1F., 1000 μm

€360.00
plus VAT &
Shipping

Product Images



Description

External measuring head for easier access to hard-to-reach measuring points Zero plate and adjustment foils included Supplied in rugged carrying case Offset-Accur: With this function, the measuring instrument can be precisely adjusted to the specific measuring range by means of a two-point calibration, thus achieving a higher precision of 1 % (or less) of the measured value Selectable units: μm , inch (mil) Auto-Power-Off SAUTER TB 2000-0.IF.: Special model for the automotive industry, precision standard 5 % of the measured value

Additional Information

No.	KP-5500
Manufacturer (Brand)	KERN
EAN	4058072150105
Transport temperature	Room temperature
autoclavable at 121°C	No
sterile	No
Manufacturer	CORE

