



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN Coating thickness gauge TE 1250-0.1FN., 1250 $\mu\text{m}$

**€520.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

External measuring head for easier access to hard-to-reach measuring points RS-232 data interface, including zero plate and adjustment foils included Supplied in rugged carrying case Offset-Accur: This feature allows the gauge to be accurately set to the specific measurement range through a two-point calibration to achieve a higher precision of 1% (or less) of the measured value Selectable units:  $\mu\text{m}$ , inch (mil) Auto Power-Off

## Additional Information

---

No.	KP-5513
Manufacturer (Brand)	KERN
EAN	4058072150235
Transport temperature	Room temperature
autoclavable at 121°C	No
sterile	No
Manufacturer	CORE

