




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Coating thickness gauge TE 1250-0.1N., 1250 μm

€450.00
plus VAT &
Shipping

Product Images



Description

External measuring head for easier access to hard-to-reach measuring points RS-232 data interface, including zero plate and adjustment foils included Supplied in rugged carrying case Offset-Accur: This feature allows the gauge to be accurately set to the specific measurement range through a two-point calibration to achieve a higher precision of 1% (or less) of the measured value Selectable units: μm , inch (mil) Auto Power-Off

Additional Information

No.	KP-5514
Manufacturer (Brand)	KERN
EAN	4058072150242
Transport temperature	Room temperature
autoclavable at 121°C	No
sterile	No
Manufacturer	CORE

