



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Coating Thickness Gauge TG 1250-0.1FN., 1250 μm

€600.00
plus VAT &
Shipping

Product Images



Description

LCD display, backlit, display of all information at a glance Offset Accur: This function allows the meter to be accurately set to the specific measurement range through a two-point calibration to achieve a higher precision of 1% (or less) of the measured value Scan mode for continuous measurements or single point measurement Mini statistics function: Displays measured value, average value, max and min value Internal data memory for up to 99 values Selectable units: μm , inch (mil) Zero plate and adjustment foils included RS-232 data interface supplied as standard Supplied in rugged carrying case External measuring head for easier access to hard-to-reach measuring points

Additional Information

No.	KP-5517
Manufacturer (Brand)	KERN
EAN	4058072150273
Transport temperature	Room temperature
autoclavable at 121°C	No
sterile	No
Manufacturer	CORE

