



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Ultrasonic material thickness gauge TB 200-0.1US-red., 0.1 mm

€345.00
plus VAT &
Shipping

Product Images



Description

External measuring head for easier access to hard-to-reach measuring points Zero plate for adjustment Including auto power-off Selectable units: mm, inch TB 200-0.IUS-RED. measures only the following materials: cast iron, aluminum, copper, brass, zinc, quartz glass, polyethylene, PVC, gray cast iron, nodular cast iron, steel Supplied in a sturdy carrying case

Additional Information

No.	KP-5504
Manufacturer (Brand)	KERN
EAN	4058072150143
Transport temperature	Room temperature
autoclavable at 121°C	No
sterile	No
Manufacturer	CORE

