



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Ultrasonic Material Thickness Gauge TD 225-0.1US., 0.1 mm

€485.00
plus VAT &
Shipping

Product Images



Description

External measuring head for easier access to hard-to-reach measuring points RS-232 data interface included Zero plate for adjustment included Selectable units: mm, inch Supplied in rugged carrying case

Additional Information

| | |
|-----------------------|------------------|
| No. | KP-5510 |
| Manufacturer (Brand) | KERN |
| EAN | 4058072150204 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| sterile | No |
| Manufacturer | CORE |

