



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN Ultrasonic Material Thickness Gauge TN 300-0.1US., 0.01 mm

**€740.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

External measuring head Data interface USB, standard (only for models with readability = 0.01 mm) Scan mode (10 measurements per second) or single measuring point selectable Internal data memory for up to 20 files (with up to 100 individual values per file) Selectable units: mm, inch Supplied in robust carrying case

## Additional Information

---

No.	KP-5535
Manufacturer (Brand)	KERN
EAN	4058072150402
Transport temperature	Room temperature
autoclavable at 121°C	No
sterile	No
Manufacturer	CORE

