



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Ultrasonic Material Thickness Gauge TN 80-0.01US., 0.01 mm

€700.00
plus VAT &
Shipping

Product Images



Description

External measuring head Data interface USB, standard (only for models with readability = 0.01 mm) Scan mode (10 measurements per second) or single measuring point selectable Internal data memory for up to 20 files (with up to 100 individual values per file) Selectable units: mm, inch Supplied in robust carrying case

Additional Information

No.	KP-5537
Manufacturer (Brand)	KERN
EAN	4058072150426
Transport temperature	Room temperature
autoclavable at 121°C	No
sterile	No
Manufacturer	CORE

