



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Ultrasonic Material Thickness Gauge TU 300-0.01US., 0.01 mm

€1,410.00
plus VAT &
Shipping

Product Images



Description

External measuring head for easier access to hard-to-reach measuring points Zero plate for adjustment Included data interface RS-232 Supplied in rugged carrying case Scan mode (10 measurements per second) or single measuring point selectable Internal data memory for up to 20 files (with up to 100 individual values per file) Limit value function: programming of max/min, in tension and compression direction, with output of an acoustic and optical signal. Ideal operating mode for efficient and error-free testing of serial parts Selectable units: mm, inch Robust metal housing

Additional Information

No.	KP-5541
Manufacturer (Brand)	KERN
EAN	4058072150464
Transport temperature	Room temperature
autoclavable at 121°C	No
sterile	No
Manufacturer	CORE

