



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN digital length measuring instrument LB 300-2., measuring range 300 mm, readability 0,01 mm

€1,250.00
plus VAT &
Shipping

Product Images



Description

Digital caliper with high measuring precision, even at high traversing speeds Easy to mount on machine tools, fixtures, test stands, etc. Zero setting, pre-draw and pre-length as well as unit change are manually adjustable Data interface RS-232 as standard Technical data Measuring range: 300 mm Readability: 0.01 mm

Additional Information

No.	KP-4675
Manufacturer (Brand)	KERN
EAN	4058072142766
Transport temperature	Room temperature
Manufacturer	CORE

