



neolab Migge GmbH

Rischerstr. 7-9

69123 Heidelberg

Deutschland

+49 (0)6221 / 8442-44

<https://www.neolab.de>

Umsatzsteuer-Identifikationsnummer:

DE 143 450 657



Julabo flow meter CORIOLIS 1 DN15 Measuring principle:

Coriolis

Price on request

plus VAT & Shipping

Product Images

Description

- Measuring principle: Coriolis - Applicable temperature control fluids: Thermal HL30, HL60, HL90, Water-Glycol, Fluorinert ® - Working temperature range: -40°C ... +150°C - Measuring accuracy: 1.10 ... 0.10 % - Flow rate: 1 ... 100 l/min 1) - Liquid connections: M30x1.5 external 2) - Pipe diameter: DN15 - Signal output: 4 ... 20 mA - Power supply: 100-230V/50-60Hz - Power supply and 2.5m analog signal cable for connection to PRESTO included. - Mounted in a rack (height approx. 500mm) 1) Possibly electronically controlled medium bypass for small flow rates necessary 2) Possibly adapter necessary

Additional Information

No.	JN-0638
Manufacturer (Brand)	Julabo
EAN	4058072846411
Transport temperature	Room temperature
Net weight	16.80000 kg
Country of origin	Switzerland
Customs tariff number	90268020

