



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Julabo flow meter CORIOLIS 1 DN25 Measuring principle: Coriolis

**Price on
request
plus VAT &
Shipping**

Product Images

Description

- Measuring principle: Coriolis - Applicable temperature control liquids: Thermal HL30, HL60, HL90, Water-Glycol, Fluorinert® - Working temperature range: -40°C ... +150°C - Measuring accuracy: 3.01 ... 0.10 % - Flow rate: 1...300 l/min 1) - Liquid connections: M30x1.5 external 2) - Pipe diameter: DN25 - Signal output: 4...20 mA - Power supply: 100-230V/50-60Hz - Power supply and 2.5m analog signal cable for connection to PRESTO included. - Mounted in a rack (height approx. 500mm) 1) Possibly electronically controlled medium bypass for small flow rates necessary 2) Possibly adapter necessary

Additional Information

No.	JN-0639
Manufacturer (Brand)	Julabo
EAN	4058072846428
Transport temperature	Room temperature
Net weight	18.50000 kg
Country of origin	Switzerland
Customs tariff number	90268020

