




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Linen tester plastic, 30 x 30 mm

€27.00
plus VAT &
Shipping

Product Images



Description

Precision linen tester with lateral scale with mm graduation.

Housing made of black plastic.

PXM® light lens, coated.

Cut-out: 30 x 30 mm

Additional Information

No.	2-2440
Manufacturer (Brand)	neoLab
EAN	4058072061012
Transport temperature	Room temperature
Color	Black
Wide	53.5 mm
Form	rectangular
Height	53 mm
Length	40 mm
Depth	40 mm
DM Lens	26 mm
Material lens	Polymethylmetacrylate (PMMA. PXM) (acrylic glass)
Ground	convex
Detail	30 x 30 mm
Material magnifier	Polyamide (PA. unspecified) glass fiber reinforced
Type magnifier	Linen tester
Magnification	5-fold

