




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® radio thermo-/hygrometer with 3 wireless transmitters

€106.50
plus VAT &
Shipping

Product Images



Description

The **neoLab wireless thermo/hygrometer with 3 wireless transmitters** is used to control indoor and outdoor climate. With the three included wireless sensors, temperature and humidity can be measured and observed up to a range of 100 meters. In the laboratory, the climate in three different temperature-sensitive rooms (cold room, boiler room,...) can thus be controlled independently.

The high-resolution display with dimensions of 128 x 32 x 128 mm and a weight of 217 g is used for permanent visual control of the desired parameters. The display is perfectly suitable for standing or hanging.

- 3 x 1.5 V AA batteries not included
- Temperature range: -10 to 60 °C
- Humidity 1 to 99

Compatible with neoLab wireless temperature monitoring system with 3 wireless transmitters (2-5442) for permanent data recording on PC or Android devices.

Additional Information

No.	2-5441
Manufacturer (Brand)	neoLab
EAN	4058072156411
Transport temperature	Room temperature
Color	Black
Temp MAX	60 °C
Temp MIN	-10 °C
Wide	128 mm
Height	32 mm
Depth	128 mm
Display (reading, setting)	digital
Weight	217 g
Type interface	WLAN
Counting and measuring devices Features	with interface with limit value setting with alarm function in set including batteries/battery suitable for wall mounting
(Combined) functions	Thermometer Hygrometer
Number of parts in set	2
Energy supply	Battery Battery
Measuring unit	mm HG rel. °C °F
measured size(s)	Temperature Humidity
rF MAX	99.00 rH
rF MIN	1.00 rH

