



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® radio temperature monitoring system with 3 wireless transmitters

**€192.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### **neoLab wireless temperature monitoring system with 3 wireless transmitters**

It has never been so easy to professionally monitor and document temperature values. The neoLab wireless temperature monitoring system with 3 wireless temperature transmitters and waterproof cable probe is the ideal solution for professional monitoring and documentation via internet browser.

Optimal for the control of indoor or outdoor temperature and simultaneous monitoring of cooling and air conditioning devices or liquids via cable probes. Samples and devices that need to be cooled permanently are conveniently under control via a PC.

- Continuous recording of measured values for seamless monitoring and documentation
- Temperature range: -30 to +60 °C
- Cable length 140 cm
- suitable for PC, Apple, Linux, Android and iOS devices

## Additional Information

No.	2-5442
Manufacturer (Brand)	neoLab
EAN	4058072156428
Transport temperature	Room temperature
Color	White
Temp MAX	60 °C
Temp MIN	-30 °C
Wide	41 mm
Height	35 mm
Depth	129 mm
Display (reading, setting)	digital
Weight	66 g
Type interface	WLAN
Counting and measuring devices Features	with interface with limit value setting with alarm function in set including batteries/battery suitable for wall mounting
(Combined) functions	Thermometer Data logger
Number of parts in set	4
Data transmission	WLAN radio
Energy supply	Battery Battery
Spare part	2-3023
Measuring unit	°C °F
measured size(s)	Temperature

