



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Moisture analyzer 0.001 g : 160 g

€1,490.00
plus VAT &
Shipping

Product Images



Description

- Backlit LCD display, digit height 17 mm
- Instantaneous moisture content in %
- Unit of result display
- Drying process active
- Halogen quartz lamp 400 W
- Scale-internal memory for automatic sequence of 5 complete drying processes
- The last measured value remains on the display until it is overwritten by a new measurement
- 10 sample pans included
- Protective work cover included
- Application manual: On the internet you will find a practical application manual for each KERN moisture analyzer with many examples, field reports, settings and tips
- Tip: Faster, more precise, more comparable results - this is how you get the most out of your KERN moisture analyzer!

Additional Information

No.	KP-4010
Manufacturer (Brand)	KERN
EAN	4058072136871
Transport temperature	Room temperature
Readability	0,001 g
Manufacturer	CORE
Material weighing plate	Stainless steel
Sample weight MAX	160 g
Weighing range MAX	160 g

