



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Moisture analyzer 1 mg : 60 g

€1,650.00
plus VAT &
Shipping

Product Images



Description

- Backlit graphic LCD display with digit height 15 mm
- Drying process active
- Drying time so far
- Current temperature
- Unit of result display
- Current moisture content in %
- Active heating profile
- Halogen quartz glass heater 400 W
- Viewing window on the sample, very helpful for the initial adjustment
- Balance internal memory for automatic running of 10 drying programs and 100 drying processes performed
- The last measured value remains on the display until it is overwritten by a new measurement
- Sample designation for up to 99 samples, 2-digit, freely programmable, is printed in the measurement report
- Date and time display as standard
- 10 sample dishes included
- Protective work cover included
- Application manual: On the internet you will find a practical application manual for each KERN moisture analyzer with many examples, field reports, settings and tips
- Tip: Faster, more precise, more comparable results – this is how you get the most out of your KERN moisture analyzer!

Additional Information

No.	KP-7291
Manufacturer (Brand)	KERN
VGKL number	105835603
old neoLab article no.	5-3071
EAN	4058072284091
Transport temperature	Room temperature
Manufacturer	CORE

