



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Climate Logger Smarthy

€46.50
plus VAT &
Shipping

Product Images



Description

Turn your smartphone into a mobile climate measuring device. The neoLab climate logger Smarthy measures temperature and humidity. Read out the data using a free smartphone app, no battery required. Climate logger to plug into the headphone jack High and low values Graphical history display Freely adjustable alarm values Dew point calculation Technical data: Temperature: -20 to +60°C Humidity: 20 to 95%rh Dimensions: 18 x 9 x 54 mm Weight: 7 g

Additional Information

No.	2-5685
Manufacturer (Brand)	neoLab
EAN	4058072069933
Transport temperature	Room temperature
Color	White
Temp MAX	60 °C
Temp MIN	-20 °C
Wide	18 mm
Height	9 mm
Depth	54 mm
Display (reading, setting)	digital
Weight	7 g
Counting and measuring devices Features (Combined) functions	with integrated data memory Thermometer Hygrometer
Resolution rH %	1.00 %
Resolution °C	0.1 °C
Energy supply	without power supply
Measuring unit	°C °F rH %
Measuring accuracy rH %	1.00 rH
Measuring accuracy °C	0.1 °C
Calibration type	automatic
measured size(s)	Temperature Humidity
rF MAX	95.00 rH
rF MIN	20.00 rH

